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Cat: G1285

Size: 6×50mL/6×100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent	6×50mL	6×100mL	Storage
Reagent (A):Alcian Blue Staining Solution	50mL	100mL	2-8°C, avoid light
Reagent (B):PAS Oxidant	50mL	100mL	2-8°C, avoid light
Reagent (C):Schiff Reagent	50mL	100mL	2-8°C, avoid light
Reagent (D):Hematoxylin Solution	50mL	100mL	RT, avoid light
Reagent (E):Acidic Differentiation Solution	50mL	100mL	RT
Reagent (F):Scott Bluing Solution	50mL	100mL	RT

Introduction

Glycogen staining is one of the conventional staining methods in pathology. McManus first used PAS technology to display mucin in 1946. This method is often used to display glycogen and other polysaccharides. The staining solution can not only display glycogen, but also show neutral mucilaginous substances and some acidic substances.

The combination of alcian blue and PAS technology can identify the neutral mucin and acid mucin in the same tissue section. This technique is also commonly used to detect mucins. The staining result of this technique is negative, which can clearly conclude that the substance is not mucin. In most methods, sections are first stained with standard alcian blue (pH 2.5) and then follow PAS technology. Alcian blue can dye salivary mucin, hiomucin and proteoglycan blue. PAS technology can dye neutral mucin into deep red or red purple, and at the same time dye tissues and cells containing both neutral mucin and acid mucin into different shades of purple, which is due to the combination and reaction of alcian blue and Schiff Reagent. The above staining often occurs in goblet cells of small intestine containing neutral mucin and salivary mucin.

Alcian blue is a kind of copper titanium cyan dye. This kind of cationic dye combines with acid group, that is, alcian blue forms insoluble complex with anion group such as carboxyl group and sulfate group in the tissue. The salt bond with positive charge in the molecule combines with the acid group with negative charge in the acidic mucopolysaccharide substance to form an insoluble complex, which is blue, and then combine with PAS for compound dyeing, so as to display three different mucilage substance components

Self Provided Materials

Distilled water, Series of ethanol

Protocol (for reference only)

1. Dewax to distilled water, then rinse in distilled water for 2min.
2. Stain with Alcian Blue Staining Solution for 10-20min.
3. Rinse in distilled water three times for each time 1-2min.
4. Treat by Oxidant for 5min. Rinse in tap water and rinse in distilled water twice.
5. Soak in Schiff Reagent and stain for 10-20min. Then wash with running water for 10min.
6. Stain with Hematoxylin Solution for 1-2min, then wash with water.
7. Differentiate by Acidic Differentiation Solution for 2-5s, then wash with water.
8. Blue with Scott Bluing Solution for 3 min and wash with water for 3 min.
9. Dehydrate by series of ethanol, transparent by xylene and seal with resinene.

Result

Glycogen, neutral mucins and various glycoproteins	Purplish Red
Acidic mucins (sulfated and carboxylated)	Blue
Proteoglycans, hyaluronic acid	Blue
Note: the cells or tissues that containing neutral mucins and acidic mucins can be dyed different shades of blue purple to purple.	

Note

1. Section dewaxing should be as clean as possible, otherwise it will affect the dyeing effect.
2. The oxidation time of Oxidant should not be too long, and the best temperature is 18-22 °C.
3. Reagent A, B and C should be kept in 4 °C airtight storage, and avoid too much sunlight and air during use. Before use, it is better to take it out in advance and restore room temperature, and use it in dark.
4. Acidic Differentiation Solution should be replaced frequently, and the differentiation time should be determined according to the thickness of section, the type of tissue and the old and new of Acid Differentiation Solution. In addition, the washing time of tap water after differentiation should be enough.
5. The time of action in Oxidant and Schiff Reagent is very important, which depends on the thickness of section and the type of tissue.
6. When re-dyeing with Hematoxylin Solution, light dyeing is required to avoid affecting the observation of positive substances. The purpose is to prevent the cytoplasm or mucin staining from covering up the color of alcian blue.
7. The results show that the dyeing sequence of Alcian Blue-PAS Stain can affect the final results. When PAS Staining is carried out before Alcian Blue Staining, neutral mucin and glycogen can be dyed purple. In contrast, Alcian Blue can dye these substances to the expected purplish red color before PAS Staining.
8. For your safety and health, please wear experimental clothes and disposable gloves.