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Cat: G1510

Size: 100mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Neutral fat stains often use Sudan II, Sudan III, Sudan IV, Sudan black B, oil red O and so on. Dyestuff is often used in traditional methods. Recently, azo dye oil red O is more suitable for dyeing fat. Oil red O is a strong lipid solvent and dye. It is easy to bind to triglycerides in droplet shape, but has a slightly poor binding to phospholipids. The dyeing principle is generally regarded as a physical miscibility or adsorption, and the fat is dyed by miscibility. The solubility of dyes in frozen sections is higher than that in the original solvent, so when dyeing, the dyes are transferred from organic solvents to lipids and the fats are dyed.

Sudan III Stain Solution is mainly used to show fatty degeneration of tissues and organs and abnormal lipid-like sedation. It often occurs in fatty degeneration of liver, kidney, heart and other parenchymal organs that most of the neutral fat droplets appear in cells. It can identify and diagnose tumors and their properties in adipose tissue. The sample should not fixed with ethanol (10% formalin if required) or paraffin. Frozen sections or carbon wax sections were needed. The positive staining results of fat were orange to red, but the specific color depended on the concentration of lipid.

Self Provided Materials

1. Slide, Microscope, Distilled Water, Glycerol Gelatin
2. 70% Ethanol, Differentiation Solution, Mayer Hematoxylin Staining Solution

Protocol (for reference only)

1. Prepare low-temperature sections of fresh tissue, generally at -20°C to -25°C. If the sample is lipoma, it should be adjusted to -30 °C.
2. Slice the frozen section at 6-15 micron (6-8 micron is preferable) and attach to the slide.
3. Wash with distilled water.
4. Re-dyeing with Mayer Hematoxylin Staining Solution for 1 min.
5. After tap water washing, differentiate by acid alcohol for several seconds, and wash with tap water till the nucleus is blue.
6. Rinse with distilled water and soak slightly with 70% ethanol.
7. Stain with Sudan III Solution for 30min.
8. Differentiate with 70% ethanol for several seconds and wash in tap water.
9. Use filter paper to absorb the slices and surrounding water and let them dry slightly.
10. Glycerol gelatin sealing.

Result

Neutral fat	Orange Red
Nucleus	Light Blue

Note

1. Samples should not be fixed with ethanol or embedded with paraffin, and frozen sections should be used. It is recommended to use 10% neutral formalin or 10% formaldehyde-calcium solution to fix.
2. Precipitation of dyes must be prevented during dyeing. Therefore, slices should be covered when they are put into dye solution, so as not to contact with flowing air to avoid precipitation when the solution volatilizes.
3. Frozen sections are easier to be colored, and over-dyeing should be avoided when re-dyeing with Mayer hematoxylin.
4. Sudan dyes fade easily and should be kept in airtight condition.
5. Samples sealed with glycerol gelatin did not last long. If long-term preservation is required, the edge of the boundary between the cover slide and the slide can be sealed with neutral gum.
6. For your safety and health, please wear lab clothes and disposable gloves.