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- [1] Yu Jiang,Dantian Zhu,Wenfeng Liu,et al. Hedgehog pathway inhibition causes primary follicle atresia and decreases female germLine stem cell proliferation capacity or stemness. Stem Cell Research & Therapy. July 2019. (IF 4.627)

## Alkaline Phosphatase Stain Kit (Kaplow's/Azo Coupling Method)

**Cat:** G1480

**Size:** 4×2mL/4×10mL/4×20mL

**Storage:** -20°C, avoid light, valid for 6 months.

### Kit Components

Reagent		4×2mL	4×10mL	4×20mL	Storage
Reagent A : ALP Fixative		2mL	10mL	20mL	RT, avoid light
Reagent B :ALP Incubation Solution	B1:AS-BI Solution	1mL	5mL	10mL	-20°C, avoid light
	B2:FBB Solution	1mL	5mL	10mL	-20°C, avoid light
Mix B1:B2 in 1:1 ratio as ALP Incubation Solution.It's ready to use.					
Reagent C :Nuclear Fast Red Solution		2mL	10mL	20mL	2-8°C, avoid light
Reagent D :Methyl Green Solution		2mL	10mL	20mL	RT, avoid light

### Introduction

Alkaline phosphatase (ALP or AKP) is widely distributed in mammalian tissues. The optimum pH for its activity is 9.2-9.8. This enzyme mainly exists in active substance exchange sites (cell membranes):Such as the brush-like margin of intestinal epithelium and proximal convoluted tubule of kidney, the stationary cilia of epididymis epithelium, the capillary cholangium of liver, and the endothelium of arterioles and capillary arteries. It is also found in the endoplasmic reticulum, Golgi complex, swallowing vesicles, lysosomes of intestinal epithelium, neutral granules of neutrophils and smooth muscle cell membranes.

Alkaline Phosphatase Staining Kit (Kaplow Method) is not a metal precipitation method to display alkaline phosphatase activity, but an azo coupling method (also known as simultaneous coupling method).The principle is that under the alkaline condition of pH 9.2-9.8, intracellular alkaline phosphatase can hydrolyze AS-BI phosphate and release phosphoric acid and naphthol, which can form colored products with coupling diazo salts and locate in the cytoplasm.The dye can be used for blood, bone marrow or cell smears, frozen and paraffin sections.The alkaline phosphatase positive site is blue and located in the cytoplasm. This method is more reliable than metal salt precipitation method.

### Self Provided Materials

Glass slide, Wet box, Optical microscope

### Protocol(for reference only)

*Reagent C :Nuclear Fast Red Solution may produce suspended solids or a small amount of precipitation due to flocculation. It is recommended to take supernatant or boil water bath for 5-10min and then air it to 30-40 °C. (see Note 2)*

#### For Smear or Section Stain

1. Fix blood, bone marrow or cell smears, frozen sections and dewaxed paraffin sections in ALP Fixative for 3min and wash with distilled water.
2. Add the ALP Incubation Solution to sections and place in a wet box. Incubate in dark for 15-20min and wash by distilled water.
3. Re-dyeing with Nuclear Fast Red Solution or Methyl Green Solution for 3-5min.
4. View the sections under microscope after washing or glycerol gelatin sealing.

#### For Adherent Cells Stain

1. Take the sample cultured in 6-well plates or other containers, remove the cell culture medium and rinse the cells completely with PBS.
2. Fix by ALP Fixative for 3mins or 4% PFA for 10-15min, then wash with PBS.
3. Add the ALP Incubation Solution to sections and place in a wet box.Incubate in dark for 15- 20min and wash with PBS.
4. Re-dyeing with Nuclear Fast Red Solution or Methyl Green Solution for 3-5min.
5. Wash with PBS and view the sections under microscope.

### Result

Positive site of ALP	Blue
Nucleus	Red or Green

**Criteria for blood and bone marrow smears judging**

Judge according to 100 neutrophil-positive granules and report the level in 0-4

Positive level	Dyeing Characteristics
0	No Granule
1	Slightly Granule
2	Medium-Grade Granule
3	Majority Granule
4	Full Granule

**Note**

1. Blood or bone marrow smears or other samples should be fresh. Cut sections in appropriate thickness and fix it in time to avoid reducing activity of the enzyme.
2. Reagent C :Nuclear Fast Red Solution is a colloidal solution, which is stored at low temperature (lower than 25 °C) or stored for a long time. Suspended solids or a small amount of precipitation are generated due to flocculation, which is a normal phenomenon and generally does not affect the use. If the colloid solution is evenly dispersed in the boiling bath, tighten the bottle cap for 5-10min to recover the turbid solution.
3. In the process of cultured cells staining, slightly wash and stain to avoid damage or loss of cells.
4. ALP Incubation Solution is easy to lose effect. Use it as soon as possible.
5. Re-dyeing by Nuclear Fast Red Solution or Methyl Green Solution.
6. There should be a positive control for each sample.

**Reference**

- [1] Yu Jiang,Dantian Zhu,Wenfeng Liu,et al. Hedgehog pathway inhibition causes primary follicle atresia and decreases female germLine stem cell proliferation capacity or stemness. Stem Cell Research & Therapy. July 2019. (IF 4.627)