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Plant Cork Cell Wall Stain Solution (Sudan Method)

Cat: G1525

Size: 50mL

Storage: RT, avoid light, valid for 1 year.

Introduction

As the defense layer of plant tissue, plant cell wall is composed of polysaccharides, proteins, lignin, lipid compounds and minerals. It provides mechanical support for tissue and plays an important role in material transportation, cell recognition and cell growth. The formation of cell wall corking is due to the increase of fatty compound cork, which makes the cell wall difficult to breathe and drain water. Finally, the Protoplast in the cell disappears and is filled by fatty compounds.

Plant Cork Cell Wall Stain Solution (Sudan Method) can be used to observe the staining of keratinized cell wall or suberized cell wall of bare handed sections, frozen sections and paraffin sections treated with plant tissue fixing solution. Keratinized cell wall or suberized cell wall has orange color reaction when treating with the solution.

Protocol (for reference only)

For frozen sections

1. **Sampling:** take the roots, stems and leaves of fresh plant tissues for cross cutting. For larger leaves, the size shall not exceed 1cm × 1cm is appropriate.
2. **Fixation:** immerse the tissue in OCT embedding agent, and then carry the embedding agent for negative pressure treatment for about 20min until the tissue is completely soaked. After quick freezing with liquid nitrogen, store at - 80 °C.
3. **Section:** the thickness of frozen section should be 10um. After drying, slightly fix tissues with precooled plant tissue fixing solution or acetone and stored at - 20 °C.
4. **Dyeing:** after rewarming, firstly treat the slices with 70% ethanol for 1-2min, soak Plant Cork Cell Wall Stain Solution (Sudan Method) for 5-10min, quickly rinse the differentiation with 70% ethanol for 2-3s, and then rinse with distilled water until there is no dye. (see note 4)
5. **Microscopic examination:** slice with water, seal with glycerin gelatin and observe under microscope.

For paraffin sections

1. **Sampling:** take the roots, stems and leaves of fresh plant tissues for cross cutting. For larger leaves, the size shall not exceed 1cm × 1cm is appropriate.
2. **Fixation:** immerse the tissue in AAF fixative for more than 4h, and then conduct negative pressure treatment with the fixative for about 20min until the tissue sinks and no gas is discharged. Transfer to 70% ethanol for equilibrium for 1h, fully elute the stationary solution, and then carry out normal dehydration, transparency and embedding treatment. (see Note 2)
3. **Production:** the thickness of paraffin section should be 3-5um, and it should be dewaxed to water normally.
4. **Dyeing:** firstly treat the slices with 70% ethanol for 1-2min, soak Plant Cork Cell Wall Stain Solution (Sudan Method) for 5-10min, quickly wash with 70% ethanol for 2-3s, and then rinse with distilled water until there is no dye. (see note 4)
5. **Microscopic examination:** slice with water, seal with glycerin gelatin and observe under microscope.

Result

Suberized Cell Wall	Orange Red
Background	Colorless or Light Yellow

Note

1. For small tissue samples (stems with diameter less than 1.5mm or leaves with thickness less than 1mm), negative pressure treatment can also be avoided.
2. If the tissue cannot be treated with negative pressure in time, can immerse tissue in AAF fixative for fixation overnight and then treat with negative pressure. Conduct the process with step 2.
3. If you want to speed up the dyeing, you can first preheat the dye at 60 °C, and the heating can promote the combination of the dye and the cork in the cell wall.
4. The reagent is volatile, and it is best to use immersion dyeing to prevent the dye from depositing on the tissue and section, which is not easy to wash off.
5. The dye can be reused for 4-5 times. When not in use, close the cover in time to prevent reagent volatilization.
6. For your safety and health, please wear experimental clothes and disposable gloves.