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Cat: G3030

Size: 4×50mL

Storage: 2-8℃, avoid light, valid for 6 months.

Kit Components

Reagent		4×50mL	Storage
Reagent (A): Gomori Silver Hexamine Solution	A1:GSH Solution A	25mL	2-8℃, avoid light
	A2:GSH Solution B	25mL	RT
Before use, mix A1 with A2 as the ratio of 1:1 to form Reagent (A), which is not suitable to make in advance.			
Reagent (B): Gold Chloride Solution		50mL	2-8℃, avoid light
Reagent (C): Hypo Solution		50mL	RT
Reagent (D): Eosin Staining Solution		50mL	RT, avoid light
Reagent (E): Gomori Control Solution		10mL	RT

Introduction

The breakdown product of purine (nucleotide) metabolism in human body is urate. Most of uric acid is not completely eliminated by the kidney, sometimes in the form of sodium urate in the blood circulation. The content of sodium urate in the blood of patients with gout is high, and the excessive sodium urate is prone to the deposition of urate, which leads to gout stone, synovitis, arthritis, kidney stone, etc. Under the microscope, most of urate are in shape of acicular crystals that parallel to each other.

Urate Stain Kit(GMS Method) is used to stain urate and make the crystal of urate black. It is often used to assist in the diagnosis of gout nodules caused by urate. Since urate is easily soluble in water but not in ethanol. When fixing, use ethanol instead of formalin or other fixatives.

Protocol(for reference only)

1. Tissue fixation: fix in absolute ethanol for 16 h or overnight. Then soak in absolute ethanol for 3 times, 30 min each time.
2. Soak in xylene twice for 20min each time, and conventionally embed in wax.
3. Cut section in 5μm thick, dewax in xylene to absolute ethanol.
4. Prepare Gomori Silver Hexamine Solution in advance and preheat it to 58-60 °C. Immerse the section in the pre heated Gomori Silver Hexamine Solution (capped) and incubate at 58-60 °C for 30min. If there is urate, the section is black. Slightly wash with distilled water.
5. Treat in Gold Chloride Solution for 1min. Slightly wash with tap water.
6. Treat in Hypo Solution for 5min. Wash with tap water for 5min.
7. Slightly stain with Eosin Staining Solution for 30s. Slightly wash with tap water.
8. Conventionally dehydrate and transparent, then seal with resinene.

Result

Urate crystal	Black
Background	Light Red

Negative control:

Take continuous sections and dewax. Incubate in Gomori Control Solution for 5min, wash twice with absolute ethanol, and then put into Gomori Silver Hexamine Solution. The other steps are the same as before, calcium salt is negative.

Note

1. Urate is easy to dissolve in water, so must use absolute ethanol for tissue fixation. Avoid the section contacting with water before entering Gomori Silver Hexamine Solution.
2. If replace thermostat with water bath, the temperature should be properly adjusted to 48-50 °C, otherwise the section will turn black.
3. If there is calcium salt in the tissue, it is easy to appear false-positive. It should be distinguished from needle like urate.
4. For your safety and health, please wear experimental clothes and disposable gloves.