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规格: 100mL/500mL

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Cat: G3668

Size: 100mL/500mL

Storage: RT, avoid light, valid for at least 1 year.

Introduction

Toluidine blue O is a kind of commonly used synthetic dye, belonging to quinone imine dyes, which is basic dye. The cations in toluidine blue have a dyeing effect. The acidic substances of tissue cells are combined with the cations in them to dye the nucleus to make it blue. The cytoplasm of mast cells contains heparin, histamine and other heterochromatic substances which encounter toluidine blue. It can be used to screen condyloma acuminatum and detect mast cells.

Protocol (for reference only)

For Tissue Slice Staining

1. Add xylene into paraffin section twice, 15min each time.
2. 1 min for each series of ethanol. Wash with tap water for 2 min.
3. Dye with Mast Cells Stain Solution (Toluidine blue Method) for 10-15mins.
4. Wash with tap water for 2 min, and drain with filter paper.
5. Differentiate with 95% ethanol till the mast cell is purple-Blue and the others are light blue.
6. Dehydration by 95% ethanol once for 1min and absolute ethanol 2 times for 2mins each time.
7. Transparent with xylene and seal with resinene.

For Cell Smear Staining

1. Dilute Toluidine Blue O Solution with 20% ethanol solution, which is generally required to be diluted to ten times.
2. After smearing the cells, fix them in 95% ethanol for 3min immediately.
3. Dye with dropping diluted Toluidine Blue O Solution, and covering the slide for 2min.
4. Put up the slide and apply a little pressure to make the excess dye be absorbed by the paper towel.
5. No need to dry, direct microscopy.

Result

Cell Nucleus, Lymphocyte	Dark Blue
Nucleolus, Mast Cell Granules	Purplish Red
Cytoplasmic, Monocyte	Light Blue

Note

1. When using this reagent for the first time, it is recommended to take 1-2 samples for pre experiment. When there are a large number of samples, a dyeing frame and a dyeing cylinder can be used for easy operation.
2. This product is an aqueous solution of 1% concentration. Dilute into specific concentration by yourself.
3. Because the temperature has a great influence on the solubility and colouring power of dyes, if you want to dye fast, you can increase the temperature properly when the room temperature is low.
4. For your safety and health, please wear experimental clothes and disposable gloves.

Related Products

- G3660 Toluidine blue staining Solution, For Cell
G3661 Toluidine Blue O solution, 1% in PBS
G3662 Toluidine Blue O Solution, 0.5% in PBS
G3663 Toluidine Blue O Solution, 1% in Sodium Borate
G3665 Toluidine Blue O Solution, 0.5% in Sodium Borate
G3670 Mast Cells Stain Solution (Toluidine blue Method)