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Elastic Fiber Stain Kit(Verhoeff Method)

Cat: G1598

Size: 2×100mL

Storage: RT, avoid light, valid for 6 months.

Kit Components

Reagent		2×100mL	Storage
Reagent(A) : Verhoeff Working Solution	Verhoeff Solution A	60mL	RT, avoid light
	Verhoeff Solution B	24mL	RT, avoid light
	Verhoeff Solution C	24mL	RT, avoid light
Before use, mix A,B and C as the ratio of 5:2:2 to prepare Verhoeff Working Solution. It is recommended to use in 2-3h			
Reagent(B): Verhoeff Differentiation Solution		100mL	RT

Introduction

In a narrow sense, connective tissue contains three types of fibers: collagen fibers, reticular fibers, elastic fibers. Elastic fibers are found in the lungs, arteries, veins. It is highly refractive, elastic and usually thinner than collagen fibers. Elastic fibers stain well with Gomori aldehyde fuchsin, orcein, resorcin-fuchsin, and Weigert's elastic stain in histological sections.

Elastic Fiber Stain Kit(Verhoeff Method) is used to demonstrate changes in elastic tissue such as atrophy, arteriosclerotic changes, reduplication, breaks, splitting, or to determine whether blood vessels have been invaded by a tumor.

Protocol(for reference only)

1. Dewax to distilled water.
2. Soak staining in Verhoeff Working Solution for 15-30 mins. Rinse in tap water.
3. Differentiate in Verhoeff Differentiation Solution, until the elastic fibers are clear. Rinse in tap water.
4. Rinse in 95% alcohol for 2-5min to make elastic fibers more clear.
5. (optional) Re-dyeing with Eosin solution or other counterstain.
6. Dehydrate in absolute alcohol and transparent in xylene, seal with resinene.

Result

Elastic Fibers	Blue purple to purple black
Other Tissue	According to the counterstain

Note

1. Most fixatives can be used to fix the sections in this dye.
2. The degree of differentiation should be controlled under microscope to avoid over staining or light staining.
3. Warm tap water washing can enhance the degree of fiber coloring.
4. For your safety and health, please wear laboratory clothes and disposable gloves.