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**Cat:** G1282

**Size:** 5×50mL/5×100mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

### Kit Components

Reagent	5×50mL	5×100mL	Storage
Reagent (A): $\alpha$ - Diastase Solution	50mL	100mL	2-8°C
Reagent (B): PAS Oxidant	50mL	100mL	2-8°C, avoid light
Reagent (C): Schiff Reagent	50mL	100mL	2-8°C, avoid light
Reagent (D): Mayer Hematoxylin Solution	50mL	100mL	2-8°C, avoid light
Reagent (E): Acidic Differentiation Solution	50mL	100mL	RT

### Introduction

Glycogen staining is one of the conventional staining methods in pathology. McManus first used PAS technology to display mucin in 1946. This method is often used to display glycogen and other polysaccharides. The staining solution can not only display glycogen, but also show neutral mucilaginous substances and some acidic substances, as well as cartilage, pituitary, mould, fungus, pigment, amyloid substance, basement membrane, etc.

The oxidant can oxidize 1,2-glycol group in sugars and related substances to make it into dialdehyde. Aldehyde and Schiff Reagent can combine to form a purplish-red fuchsin complex. As oxidants can also oxidize other substances in cells, we should pay attention to the concentration and oxidation time of oxidants when using, so that the oxidation control can oxidize glycol group to aldehyde group without peroxide, which is a key step.

The characteristic of Periodic Acid Schiff Diastase (D-PAS) Stain Kit is that before glycogen PAS staining, it is treated by diastase. When glycogen is digested, prepare two identical sections. After dewaxing, one section is treated with appropriate buffer containing diastase, and the other one is treated only with buffer. Then the two sections are stained with PAS method, and the disappearance of staining after digestion indicated the presence of glycogen.

### Protocol (for reference only)

1. Prepare two identical sections and dewax to distilled water.
2. Treat one section into  $\alpha$  - Diastase Solution for 20min. The other one directly into distilled water for 1h without treated by  $\alpha$  - Diastase Solution.
3. Rinse the two sections in tap water respectively for 5-10min.
4. Place the sections into Oxidant and incubate for 5-8min at room temperature, generally no more than 10min. Rinse in tap water once and then rinse in distilled water twice.
5. Place the sections into Schiff Reagent and incubate in dark for 10-20min. Rinse in distilled water for 10min.
6. Stain with Mayer Hematoxylin Solution for 1-2min.
7. (optional) differentiate by Acidic Differentiation Solution for 2-5s. Rinse in tap water for 10-15min to blue.
8. Dehydrate, transparent by xylene and seal with resinene.

### Result

Glycogen, Neutral and Salivary Mucin	Red to Purple
Various Glycoprotein	Red to Purple
Nucleus	Blue
The section untreated, glycogen is bright red or red purple; The section treated by Diastase Solution, glycogen shows negative.	

### Note

1. Section dewaxing should be as clean as possible, otherwise it will affect the dyeing effect.
2. A positive control section is required to verify the enzyme activity.
3. The oxidation time of Oxidant should not be too long, and the best temperature is 18-22 °C.
4. Reagent A, B and C should be kept in 4 °C airtight storage, and avoid too much sunlight and air during use. Before use, it is better to take it out in advance and restore room temperature, and use it in dark.

Acid Differentiation Solution should be replaced frequently, and the differentiation time should be determined according to the thickness of section, the type of tissue and the old and new of Acid Differentiation Solution. In addition, the washing time of tap water after differentiation should be enough.

5. The time of action in Oxidant and Schiff Reagent is very important, which depends on the thickness of section and the type of tissue.
6. The staining time of frozen section should be as short as possible.
7. For your safety and health, please wear experimental clothes and disposable gloves.