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**Cat:** G2541

**Size:** 2×50mL

**Storage:** 2-8°C, avoid light, valid for 1 year.

### Kit Components

Reagent	2×50mL	Storage
Reagent (A): Alcian Acidic Solution	50mL	RT
Reagent (B): Alcian Staining Solution(pH1.0)	50mL	2-8°C, avoid light

### Introduction

Alcian Blue 8GX also known as Alcian Blue, is a kind of copper titanium cyan conjugate dye, which was originally used for textile fiber dyeing. This kind of cationic dye combines with acid group, that is, alcian blue forms insoluble complex with anionic group such as carboxyl group and sulfate group. Alcian Blue is composed of the central copper containing phthalocyanine ring and four isothiurea groups linked by thioether bond. The isothiurea group is moderately alkaline, making Alcian Blue with cation.

Alcian Blue Cartilage Stain Kit, pH 1.0 can distinguish the genus of mucus by different pH value of staining solution. When pH value is 1, carboxyl group(COOH) is not dyed, and sulfuric acid group(OSO<sub>3</sub>H) is dyed. When pH value is 2.5, COOH is dyed well, while OSO<sub>3</sub>H is dyed bad, which makes mucin sulfate and mucin saliva stain. This reagent is only suitable for scientific research, not for clinical diagnosis or other purposes.

### Protocols(for reference only)

1. Conventional decalcification, dewax to distilled water.
2. Soak in Alcian Acidic Solution for 3min.
3. Stain with Alcian Staining Solution(pH1.0) for 30min.
4. Rinse in running water for 5min.
5. Dehydrate in series of ethanol, transparent in xylene and seal with resinene.

### Result

Sulfated Mucin	Blue
Nucleus	Light Blue

### Note

1. The fixative adopts 10% neutral formalin.
2. This method can distinguish sulfated mucin and proteoglycan when pH value is 1.0.
3. The opened reagent shall be used up within 6 months after open, and the cap shall be tightened timely after each use to avoid volatilization or deterioration.
4. For your safety and health, please wear experimental clothes and disposable gloves.

### Reference