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## Alkaline Phosphatase-PAS Stain Kit

**Cat:** G2450

**Size:** 6×50mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

### Kit Components

Reagent	6×50mL	Storage
Reagent(A): ALP Incubation Solution	50mL	2-8°C, avoid light
Reagent(B): Co Solution	50mL	RT, avoid light
Reagent(C): ALP Vulcanizing Solution	2×1mL	RT, avoid light
Reagent(D): Oxidant	50mL	2-8°C, avoid light
Reagent(E): Schiff Reagent	50mL	2-8°C, avoid light
Reagent(F): Sulphite Solution	50mL	RT, avoid light
Reagent(G): ALP Control Solution	10mL	2-8°C, avoid light

### Introduction

Alkaline phosphatase (ALP or AKP) is a kind of phosphatase, which is widely distributed in mammalian tissues, and the optimal pH of its activity is 9.2-9.8. This enzyme is mainly found in the active sites of material exchange (cell membrane), such as the brush border of intestinal epithelium and renal proximal convoluted tubules, the static cilia of epididymal epithelium, the membrane of bile capillaries of liver, the endothelium of arterioles and capillaries, and also in the endoplasmic reticulum, golgi complex, swallowing vesicles, lysosomes of intestinal epithelium, neutrophils of neutrophils, and the membrane of smooth muscle.

In 1946, McManus first used PAS method to display mucin, which is often used to display glycogen and other polysaccharides. The oxidant can oxidize the 1,2-glycol group in sugars and related substances to make it dialdehyde. Aldehyde and Schiff reagent can combine to form a fuchsin compound to produce purplish red. Because oxidants can also oxidize other substances in cells, we should pay attention to the selection of concentration and oxidation time, so that the oxidation can be controlled not only to oxidize glycol group to aldehyde group, but also to avoid peroxidation.

ALP-PAS Stain kit can not only show the activity site of alkaline phosphatase and glycogen, but also distinguish them. Frozen section and paraffin section can be used.

### Self Provided Materials

Distilled water, Incubator or Water bath

### Protocol(for reference only)

1. For paraffin section, dewax in xylene and dehydrate in ethanol, then put into water. For frozen section, directly put into water.
2. Fix frozen section in the mixture of equal amount acetone and chloroform were for 2-5min at 4 °C.
3. Take the section into ALP Incubation Solution (take negative control section into ALP Control Solution) and place it in 37 °C incubator. The frozen section is incubated for 5-15min and paraffin section for 2-12h. Rinse with running water for 2min and then put into distilled water.
4. Add Co Solution and stain in 37 °C incubator for 5min. Rinse with running water for 5min and then put into distilled water.

5. In the above process, dilute Reagent (C) with distilled water to form ALP Vulcanizing Working Solution. It is ready to use. Take the section into ALP Vulcanizing Working Solution and incubate for 1-2min.
6. Rinse with running water for 10min and then put into distilled water.
7. Oxidize with oxidant at room temperature for 5-10min. Rinse with distilled water.
8. Place the section into Schiff Reagent and dye in the incubator for 10-15min.
9. Prepare Sulphite Working Solution according to the ratio of Sulphite Solution: distilled water is 1:10. Wash with running water for 2min and wash with distilled water for 5min .
10. For frozen section, directly seal with glycerin gelatin. For paraffin section, conventionally dehydrate in ethanol and transparent by xylene, then seal with resinene.

## Result

Alkaline Phosphatase Active Site	Black
PAS Positive Substance	Purplish Red

### Positive Control(*optional*):

1. ALP Control Solution is an incubation solution without substrate. Take the same section into ALP Control Solution, the rest follow the above steps. The result is negative.
2. (*Alternative*) Before the section enters into the ALP Incubation Solution, it can be rinsed by Iodine Solution and 5% Sodium Thiosulfate Solution for 3min respectively, and wash with water fully. Then it can be incubated and continue the rest steps. This method can be used as a negative control.

## Note

1. The staining solution is suitable for frozen section and paraffin section. Pay attention to control the treatment time of oxidant in dyeing process, otherwise it is easy to fade.
2. After the color development of alkaline phosphatase, pay more attention after oxidation by oxidant, and control the time strictly, otherwise the color will fade.
3. ALP Incubation Solution and ALP Vulcanizing Solution are easy to lose effect, so it is better to store them in small parts and use them immediately once open.
4. ALP Vulcanizing Solution has corrosive and irritating smell, so operate carefully.
5. The exposure time of frozen section at room temperature should be reduced when staining frozen section.