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规格: 100mL /500mL

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Cat: G1613

Size: 100mL /500mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Papanicolaou method is widely used in routine cytological staining. Papanicolaou stain is initially used to detect estrogen levels in vaginal epithelium and pathogens such as Candida and trichomonad in the genital tract. The combination of orange G6 and EA36 or EA50 can dye the cytoplasm into bright green, blue and pink. At present, most laboratories use finished dye solution, so we should pay attention to the best conditions of each dye solution after improvement. The final cytoplasmic staining should be transparent, and the nuclear chromatin should be easily distinguished.

At present, the improved pasteurization solution contains many kinds of ions, which has polychromatic dyeing efficiency. After staining, the cytoplasm is bright and transparent, and the structure of nuclear membrane, nucleolus and chromatin is clear. The nuclear staining solution is mainly Harris hematoxylin staining solution, and the cytoplasm staining solution is mainly EA36 staining solution and EA50 staining solution. Papanicolaou Stain is used for exfoliated cells. The nucleus is blue or black, and the cytoplasm of keratinized squamous cells is pink or orange red.

Papanicolaou EA36 Stain Solution can make the cytoplasm of different degrees of keratinization have different colors. It is more suitable than EA50 for gynecological cytology smear staining, such as screening cervical cancer and precancerous lesions.

Self Provided Materials

Fixative (like AF Fixative), Series of alcohol, Microscope, Acid alcohol differentiation solution

Protocol (for reference only)

1. For cell smear, fix with 95% alcohol for 10-15mins.
2. (optional) Rinse in 95%, 85%, 75% alcohol separately for 1min. Rinse with tap or distilled water for 1min.
3. (optional) Stain with Hematoxylin Solution for 5-10mins. Rinse with tap water for 2mins.
4. (optional) Differentiate with 1% acid alcohol differentiation solution for about 4-5s or 0.5% acid alcohol differentiation solution for about 10s. Rinse with tap water for 2mins.
5. (optional) Blue with Bluing Solution for 2mins. Rinse with tap water for 2mins.
6. (optional) Dehydrate in 75%, 85%, 95% alcohol separately for 2mins.
7. (optional) Stain with Orange G Solution for 2mins. Dehydrate in 95% alcohol for 2mins twice.
8. Stain with EA36 Solution for 2mins. Dehydrate in 95% alcohol for 2mins twice.
9. Dehydrate in absolute alcohol for 1min twice.
10. Transparent with xylene and seal with resinene.

Result

Nucleus	Blue Purple or Black
Cytoplasm of non keratinized cells	Light Blue or Light Green
Cytoplasm of keratinocytes	Pink or Orange

Note

1. The dye solution shall be filtered before use. If it is reused, the dye solution shall be changed frequently.
2. This product contains a certain proportion of alcohol components. Pay attention to capping and dyeing to prevent volatilization. At the same time, pay attention to capping in time for large packaging.
3. For your safety and health, please wear experimental clothes and disposable gloves.