产品合格证 Product Conformity Certificate

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1. Product Brief

L-VM-MINI vortex mixer is equipped with a upgraded brushless DC motor, allowing 4000rpm at full speed for vigorous and powerful mixing action. Plus, 5 seconds rapid mixing, mixing test tubes up to 50mL. This mini vortex shaker provides bigger area to contact with the mixing bottle, fits various types of containers. Easy one-handed touch operation, the vortex mixer starts working when slight pressure is applied to the center of the silicon cap. Low activation force, easy for repetitive and continuous mixing.

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3. Product Data

Model	L-VM-MINI
Operation Mode	Touch Operation
Voltage	100-240V 50/60HZ
Motor Power	5W
Material	ABS Plastic + Silicone
Mixing Speed	4000rpm (Fixed)
Max. Mixing Capacity	50mL
Motor Type	Brushless DC Motor
Type of Movement	Orbital
Orbital Diameter	20 mm
Protction class	IP42
Permissible Ambient Temperature	41°F-104°F / 5°C-40°C
Permissible Relative Humidity	80%
Weight	0.83 lbs / 380g
Size	3.9" x 3.9" x 2.8" / 10 x 10 x 7 cm

4. Application

It is mainly used for molecular hybrid blood and other samples that need to be mixed.

5. Usage

It is widely used as an experimental tool in institutions of higher learning.
environmental protection, sanitation, epidemic prevention, chemical industry, tap water and medical treatment.

6. Steps for Usage

- > Plug it in and turn the knob to the right to adjust the speed;
- > Place the centrifuge tube or test tube to the top center, click mode, touch and press lightly to run.

7. Instrument Using Environment

- > Indoor use
- > Altitude: ≤2000m
- > The working temperature range of the instrument is +5°C~+40°C
- > The applicable temperature range of the instrument is ≤80%
- > There is no vibration and airflow that affect performance around

8. Safety Precautions

- > Do not plug or unplug the power connector and toggle the power button when the hands are stained with liquid.
- > Do not plug or unplug the power plug when the instrument is powered.
- > It is strictly forbidden to maintain and scrub the instrument under power-on state.
- > It is strictly forbidden to install the instrument on a work surface that is uneven and shakes and shakes.